





ALL ABOUT DEVERRA HA FILLEERS

Hyaluronic acid (HA) is a naturally occurring substance that keeps skin moisturized and soft.

DEVERRA fillers are nonsurgical options that contain a modified form of HA.

Treatment with DEVERRA products by a licensed specialist provides long-lasting, natural-looking results.

DEVERRA is the hyaluronic acid fillers.

Each product in the DEVERRA Collection of Fillers adds volume to a different area of the face and body without surgery.

All DEVERRA fillers are smooth, crystal-clear gels that are administered by injection in your aesthetic healthcare provider's Clinic.

Commitment is to provide patients a wide range of new medical solutions with an excellent of safety and efficacy.

We promise to be the Brand with trust and faith by providing high quality products that give the best value.

Medical devices for professional use.

Only to be injected by an authorized medical professional.



Safety

Efficacy

Long lasting

High Viscosity


- The material is biodegradable and biocompatible
- It is absorbed in the body
 - It stimulates collagen and restore the lost volume.
 - Made from a non-animal origin stabilised hyaluronic acid.




Name	Figures	Product specification							Needle specification	
		HA	PHASIC (GEL TYPE)	BDDE	CONTENTS	LIDOCAINE	VISCOICITY COMPLEX	INDICATIONS	GAUGE (UTW)	QUANTITY
FINE PLUS Lidocaine	DEVERRA Cross Linked Hyaluronic Acid 1ml FINE PLUS With Lidocaine	24 mg/ml	MONO	<0.2ppm	1ml syringe	0.3%	180~199 Pa	Superficial dermis	30G	1EA
UNIVERSAL PLUS Lidocaine	DEVERRA Cross Linked Hyaluronic Acid UNIVERSAL PLUS With Lidocaine	24 mg/ml	MONO	<0.2ppm	1ml syringe	0.3%	280~299 Pa	Mid dermis, Subcutaneous	27G/30G	1EA
INTENSE PLUS Lidocaine	DEVERRA Cross Linked Hyaluronic Acid IMI INTENSE PLUS With Lidocaine	24 mg/ml	MONO	<0.2ppm	1ml syringe	0.3%	300~370 Pa	Mid/Deep dermis, Subcutaneous	27G	1EA
VOLUME PLUS Lidocaine	DEVERRA Cross Linked Hyaluronic Acid 1ml VOLUME PLUS With Lidocaine	24 mg/ml	MONO	<0.2ppm	1ml syringe	0.3%	400~450 Pa	Deep Subcutaneous	25G/27G	1EA



- Improving wrinkle
- Volumizing, Contouring
- Corrective treatment for asymmetry
- Improving skin quality





BEFORE AFTER



- FINE PLUS
- UNIVERSAL PLUS
- **INTENSE PLUS**
- **VOLUME PLUS**





Treatment of fine lines wrinkles
Around and under eyes
Smoothing of the forehead wrinkles
Crease treatment
Wrinkles of periorbital area

Composition:


Cross-Linked Hyaluronic Acid 24mg/ml Lidocaine HCI 0.3% 1ML prefilled syringes


Results last:

12 - 18 Months

Packaging:

1 ML prefilled syringes2 Sterile needles (30G)IFU







Best for add volume to lips
Improve lips definition and shape
Treatment of medium wrinkles

Composition:


Cross-Linked Hyaluronic Acid 24mg/ml Lidocaine HCI 0.3% 1 ML prefilled syringes)


Results last:

12 - 18 Months

Packaging:

1 ML prefilled syringes2 Sterile needles (27G, 30G)IFU







Treatment of deep wrinkles
Add volume to lips definition shape
Nose contouring

Composition:


Cross-Linked Hyaluronic Acid 24mg/ml Lidocaine HCI 0.3% 1 ML prefilled syringes)

Results last:

12 - 18 Months

Packaging:

1 ML prefilled syringes2 Sterile needles (27G, 30G)IFU







Add extra volume to cheek, chin, jawline, face contouring

Composition:


Cross-Linked Hyaluronic Acid 24mg/ml Lidocaine HCI 0.3% 1 ML prefilled syringes)


Results last:

12 - 18 Months

Packaging:

1 ML prefilled syringes2 Sterile needles (27G)IFU





DEWERRA

Cross Linked Hyaluronic Acid





2F/205, 206 Alpha City 1-Ro, Suseong-Gu, Republic of Korea (42250) Web: WWW.DEVERRA.CO.KR TEL. +82-53-217-7302 FAX. +82-508-941-5448 Made in Korea

