

DEVERRA

Cross Linked Hyaluronic Acid

I AM DEVERRA





DEVERRA

Cross Linked Hyaluronic Acid

ALL ABOUT **DEVERRA** HA FILLERS

Hyaluronic acid (HA) is a naturally occurring substance that keeps skin moisturized and soft.

DEVERRA fillers are nonsurgical options that contain a modified form of HA.

Treatment with DEVERRA products by a licensed specialist provides long-lasting, natural-looking results.

DEVERRA is the hyaluronic acid fillers.

Each product in the DEVERRA Collection of Fillers adds volume to a different area of the face and body without surgery.

All DEVERRA fillers are smooth, crystal-clear gels that are administered by injection in your aesthetic healthcare provider's Clinic.

Commitment is to provide patients a wide range of new medical solutions with an excellent of safety and efficacy.

We promise to be the Brand with trust and faith by providing high quality products that give the best value.

Medical devices for professional use.

Only to be injected by an authorized medical professional.

DEVERRA

Cross Linked Hyaluronic Acid

Safety

Efficacy

Long lasting

High Viscosity


- The material is biodegradable and biocompatible
- It is absorbed in the body
- It stimulates collagen and restore the lost volume.
- Made from a non-animal origin stabilised hyaluronic acid.



Name	Figures	Product specification							Needle specification	
		HA	PHASIC (GEL TYPE)	BDDE	CONTENTS	LIDOCAINE	VISCOICITY COMPLEX	INDICATIONS	GAUGE (UTW)	QUANTITY
FINE PLUS Lidocaine		24 mg/ml	MONO	<0.2ppm	1ml*2syringe	0.3%	180~199 Pa	Superficial dermis	30G	2EA
UNIVERSAL PLUS Lidocaine		24 mg/ml	MONO	<0.2ppm	1ml*2syringe	0.3%	280~299 Pa	Mid dermis, Subcutaneous	27G/30G	2EA
INTENSE PLUS Lidocaine		24 mg/ml	MONO	<0.2ppm	1ml*2syringe	0.3%	300~370 Pa	Mid/Deep dermis, Subcutaneous	27G	2EA
VOLUME PLUS Lidocaine		24 mg/ml	MONO	<0.2ppm	1ml*2syringe	0.3%	400~450 Pa	Deep Subcutaneous	25G/27G	2EA


★ THE PURPOSE OF DEVERRA FILLERS

- Improving wrinkle
- Volumizing, Contouring
- Corrective treatment for asymmetry
- Improving skin quality



BEFORE


Fillers injected into the skin to help fill in facial wrinkles



AFTER

★ THE INJECTION AREA

- FINE PLUS
- UNIVERSAL PLUS
- INTENSE PLUS
- VOLUME PLUS



DEVERRA

FINE PLUS With Lidocaine

Cross Linked Hyaluronic Acid

Used For :

Treatment of fine lines wrinkles
Around and under eyes
Smoothing of the forehead wrinkles
Crease treatment
Wrinkles of periorbital area

Composition :

Cross-Linked Hyaluronic Acid 24mg/ml
Lidocaine HCl 0.3%
2ML(1 ML x 2 prefilled syringes)

Results last :

12 - 18 Months

Packaging :

1 ML × 2 prefilled syringes
4 Sterile needles (30G)
Brochure



DEVERRA

Cross Linked Hyaluronic Acid

UNIVERSAL PLUS With Lidocaine

Used For :

Best for add volume to lips
Improve lips definition and shape
Treatment of medium wrinkles

Composition :

Cross-Linked Hyaluronic Acid 24mg/ml
Lidocaine HCl 0.3%
2ML(1 ML x 2 prefilled syringes)

Results last :

12 - 18 Months

Packaging :

1 ML × 2 prefilled syringes
4 Sterile needles (27G, 30G)
Brochure



DEVERRA

INTENSE PLUS With Lidocaine

Cross Linked Hyaluronic Acid

Used For :

Treatment of deep wrinkles

Add volume to lips definition shape

Nose contouring

Composition :

Cross-Linked Hyaluronic Acid 24mg/ml

Lidocaine HCl 0.3%

2ML(1 ML x 2 prefilled syringes)

Results last :

12 - 18 Months

Packaging :

1 ML × 2 prefilled syringes

4 Sterile needles (27G)

Brochure



DEVERRA

VOLUME PLUS With Lidocaine

Cross Linked Hyaluronic Acid

Used For :

Add extra volume to cheek, chin, jawline, face contouring

Composition :

Cross-Linked Hyaluronic Acid 24mg/ml

Lidocaine HCl 0.3%

2ML(1 ML x 2 prefilled syringes)

Results last :

12 - 18 Months

Packaging :

1 ML × 2 prefilled syringes

4 Sterile needles (25G, 27G)

Brochure



DE\\VERRA

Cross Linked Hyaluronic Acid



NEWZONE

2F/205, 206 Alpha City 1-Ro, Suseong-Gu, Republic of Korea (42250)
Web: WWW.DEVERRA.CO.KR
TEL. +82-53-217-7302 FAX. +82-508-941-5448

Made in Korea

